

A close-up photograph of a woman's shoulder and upper chest. She has wet skin with water droplets. A white, circular dressing is applied to her shoulder, covering a small wound. The background is a soft, out-of-focus light blue.

Clinically approved!

Can be left in place
for up to 14 days

Leukoplast®
Leukomed® Control
Clearly controlling the risk of infection

Leukomed®
by **Leukoplast®**

Leukomed® Control

Leukomed® Control is indicated for dry and low exuding acute wounds such as surgical incision, lacerations and abrasions.



Clear wound visibility



Early detection of infection



Secure bacteria-proof barrier

Transparent, conformable film

- Enables assessment of visible signs of infection
- Reduces the need for unnecessary dressing changes

Bacteria and virus waterproof barrier

- Protects against water, bacteria, viruses and contaminants
- Helps to reduce the risk of SSI

Skin friendly adhesive

- Wide adhesive edge for high durability
- Secure yet easy to remove

Absorbent hydropolymer gel

- Controlled exudate absorption supports moist wound healing
- Pad design aids conformability



Leukomed® Control

Ref-no.	Size	Pad Size	Items per unit
73230-00	5 cm x 7 cm	2 cm x 4 cm	10
73230-01	7 cm x 10 cm	4 cm x 7 cm	10
73230-02	8 cm x 15 cm	4.5 cm x 11.5 cm	10
73230-03	10 cm x 24 cm	6 cm x 20 cm	5
73230-05	10 cm x 35 cm	6 cm x 31 cm	5

Quick mobilisation of patients

A declared aim of many caregivers is the rapid discharge of patients home or to rehab after surgery.

Leukomed® Control supports these efforts:

- Transparency offers patients the possibility to observe their wound on their own
- Flexibility and elasticity promotes everyday coping
- High quality film allows showering and washing